

## Atrazine

RaPID Assay® is a rapid field or laboratory enzyme immunoassay method for the analysis of soil and water for remediation, assessment, and industrial testing.

### Test result type

- Quantitative, semi-quantitative or qualitative depending on protocol selected

### Samples per kit

- 100 test kit (tests up to 80+ samples)

### Assay range

- 0.04 to 5.0 ppb
- LDD = 0.1 ppb

### Sample methods

- Methods exist for soil and a wide variety of commodities

### Sampling time

- 'Sample-to-Data' in approximately 60 minutes

### Specificity

- The proprietary antibody used in this assay is specific for atrazine. The cross-reactivity of the Atrazine RaPID Assay® for various triazine analogues can be expressed as least detectable dose (LDD) which is estimated at 90% B/B<sub>0</sub> or as the concentration estimated (inhibitory concentration) at 50% B/B<sub>0</sub> (IC<sub>50</sub>).

- » Rapid field testing procedure for analysis of soil and water samples
- » Quantitative data results with excellent analytical precision
- » Results available in approximately 60 minutes
- » Tests up to 50 samples at once
- » Magnetic particle immunoassay
- » Training available
- » EPA SW-846 method #4670 for drinking water



## ATRAZINE RAPID ASSAY (PPB)

Contaminant	LDD	IC <sub>50</sub>
Atrazine	0.046	0.72
Ametryn	0.053	0.39
Propazine	0.033	0.74
Prometryn	0.054	0.64
Prometon	0.056	2.22
Simazine	0.340	4.90
Terbutryn	0.090	5.50
Terbutylazine	0.310	15.5
Desethyl Atrazine	0.062	3.21
Cyanazine	1.000	>10,000
6-Hydroxy Atrazine	1.100	148
Desisopropyl Atrazine	0.800	217

### Basic Test Procedure

- Add 200 µL of calibrator or prepared sample, 250µL enzyme conjugate and 500 uL antibody coupled magnetic particles. Vortex
- Incubate for 15 minutes
- Using the RaPID magnetic separator, decant and wash
- Add 500 µL colour solution
- Incubate for 20 minutes
- Stop the reaction and read colour at 450 nm
- Quantitive results and QC parameters are calculated and printed automatically using the SmartReader Photometer and Rugged Tablet.

### Test Kit Components

- » Antibody coated magnetic particles, enzyme conjugate, colour development, stop and wash reagents for analysis of 100 test tubes
- » Standards for 0.1, 1.0, and 5.0 ppb Atrazine
- » Zero standard and 3.0 ppb control standard
- » Disposable test tubes
- » Test kit instructions

### Storage & Precautions

- » Shelf life is typically one year from date of manufacture, with specific kit expiration date information provided on product packaging
- » Reagents must be stored at 39° to 46°F (4° to 8°C) when not in use
- » Storage at ambient temperature 64° to 81°F (18° to 27°C) is acceptable for day of use
- » Kits must be brought to 64° to 81°F (18° to 27°C) before use
- » Do not expose colour solution to direct sunlight

### Required Test Materials

#### Part #

- » Atrazine RaPID Assay A00071
- » Atrazine sample diluent (100 mL) A00051

### Required Test Equipment

- » Accessory kit (purchase) 6050100A
- » Accessory kit (rental) 6997010
- » Accessory kit components:
  - » Magnetic separation rack A00004
  - » Repeater pipette A00008
  - » Adjustable volume pipette A00176
  - » SmartReader w/Tablet 6050105
  - » Vortex mixer A00014
  - » Digital timer A00015

p = purchase, r = rental