

MODERNWATER

# Multisensor2000

## On-line Trihalomethanes Monitor

The Multisensor2000 is a total Trihalomethanes (THMs) monitor system designed for potable water applications. It provides on-line, real time monitoring and alarm generation without the need for additional chemical reagents. There is no contact between the sensor and the water and virtually no moving parts, eliminating the need for frequent maintenance.

Trihalomethanes are formed as a by-product when chlorine is used to treat drinking water and they are often referred to as disinfection by-products. They result from the reaction of chlorine and/or bromine with residual organic matter present in the water being treated. THMs are potentially carcinogenic and many governments set limits on the amount permissible in drinking water.

The Multisensor2000 detects chloroform, bromoform, bromodichloromethane, and dibromochloromethane, reporting the result as total THMs in ppb ( $\mu\text{g l}^{-1}$ )

The instrument samples air from the head space in the sample tank and this non-invasive method of operation provides for reliable operation and low maintenance.

The Multisensor2000 includes the measurement instrument which is housed in a robust steel enclosure, powered from either 90 - 240V AC or 24V DC sources. Also included is a stainless steel sample tank and associated pipework and control valves. The whole system is mounted on two 12mm thick PVC boards. All that is required is the water supply, drain and power.

- No reagents
- No sensor-water contact
- Low maintenance
- Low ppb ( $\mu\text{g l}^{-1}$ ) detection limits
- Output in ppb ( $\mu\text{g l}^{-1}$ )
- Fully quantitative concentration monitoring
- Real time, online



MODERNWATER

## SPECIFICATIONS

Conforms to	UL 61010-1 / EN 61010-1:2010 EN61326-1:2006
Dimensions	Height 1460mm x width 500mm x depth 280mm
Dynamic Range	1 ppb – 1000 ppb ( $\mu\text{g l}^{-1}$ ) in water
Result Output	ppb ( $\mu\text{g l}^{-1}$ ) Total THMs
Absolute Accuracy	For measurement of 80 ppb Total THMs in water: $\pm 10\%$
Repeatability	Repeated measurement of 80 ppb Total THMs in water: $\pm 2\%$
Analysis Time	20 minutes minimum cycle time
Operating Temperature	0 - 40°C ambient air (0-30°C ambient water)
Calibration Period	Six months using a calibration kit available from Modern Water
Consumables	Every 6 months: Active carbon air filter contents / Dust filter Every 24 months: Air pump
Communications	USB, Analog 2-40 mA, Profibus, Relays (GSM modem optional)
Data Storage	SD card
Water Connections	Inlet: 1/2" BSP, 1.0 - 0.5 l/min Outlet: 1/2" BSP, 1.0 - 0.5 l/min Waste: 15mm, connected to tank through drain valve which prevents overflowing of tank

## Applications


Protection of the drinking water distribution network  
Water treatment process monitoring  
Regulatory compliance guidance  
Monitoring of swimming pools

## Process explained

The Multisensor monitors work on the principles of headspace gas equilibrium (Henry's Gas Law). Water flows through the sample tank. Equilibrium is reached in the sample tank and the gas in the headspace is sampled by an "e-nose".

The "e-nose" is an array of four proprietary sensors. Each sensor has a specific non-conductive polymer coating that selectively reacts to compounds dependent upon molecular weight, polarity and molecular diameter.

Results of THMs in the water are calculated compensating for different ambient environmental conditions.



**MODERNWATER**

tel: +44 (0)1483 696 000  
www.modernwater.com  
info@modernwater.co.uk